MANUAL DE INSTRUCCIONES

Mx Onda

IRRIGADOR DENTAL PORTÁTIL

MODELO: MX-HD2410

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MODEL: MX-HD2410

CE

PLEASE READ THIS MANUAL BEFORE OPERATION

Waste electrical products must not be disposed of with household waste.

This equipment should be taken to your local recycling centre for safe treatment.

This product complies with European Directives RoHS (2011/65/UE), on the restriction of use of certain dangerous substances in electrical and electronic appliances, transposed to the Spanish legislation by means of Real Decree 208/2005.

INTRODUCTION

MX ONDA Thank You deference has had to purchase this product. Use of this dental water jet, will help you keep in perfect health both your teeth and gums.

BASIC DENTAL CARE

Failure to clean only the surface of the teeth, and no gums or interdental space, other than cause halitosis, causes the accumulation of plaque and tartar on the teeth. Continued poor dental hygiene can lead to periodontal diseases such as gingivitis or periodontitis.

Given the importance of proper dental cleaning, the continued use of this device after every meal, it is a first step to get a perfect dental hygiene buco. With the dental irrigator can effectively clean the teeth, the interdental space, the gumline and if any, the ortodoncias.

Thanks to the pressure with which the water is projected, while cleaning the teeth produce a pleasant massage of the gums, helping improving blood circulation and preventing dental diseases.

If a slight bleeding gums occurs during the use of this product it means you could have a periodontal problem. Continue to use a couple of weeks this irrigator, bleeding gums should disappear, otherwise stop using this product and consult your usual dentist.

If the user has implemented an orthodontic appliance, it is recommended that you clean in a way, even more thoroughly, than other users.

For proper dental care buco to use this device, after each meal, you should also visit a regular basis to your regular dentist.

BEFORE CONNECTING THE DENTAL WATER JET

When using electrical appliances, they should always be followed basic safety measures to reduce the risk of fire, electric shock and / or personal injury. These measures include the following:

- Read these instructions and save them for later use carefully.
- After removing the packaging, check that the appliance is not damaged. If in doubt use it; contact the technical support staff.
- Packaging materials (plastic bags, etc.) must not be left within reach of children as they are a potential source of danger.
- If you are following some type of oral treatment or suffer from any periodontal disease, talk with your regular dentist before using this device.
- To prevent battery drain and extend the life of it, when not using the portable dental irrigator for an extended period of time, you must perform a full charge at least once a month.

SECURITY MEASURES

- This device cannot be used by children or children under 8 years or people with or without experience and knowledge reduced physical, sensory or mental capabilities, unless they have received appropriate instructions on the use of the device and be supervised by an adult or responsible person .

- This appliance is designed for household use only MX ONDA will not accept any responsibility and warranty is void in case of misuse or failure to follow instructions.
- Do not insert the device in water or any other liquid, nor leave it exposed to the weather (rain, sun, etc.).
- Do not use any accessory not recommended by the manufacturer.
- Do not use if damaged or clogged nozzles.
- Do not allow children to use this appliance without adult supervision.
- If the dental water jet is used by several people, for hygiene reasons it is recommended that each person use your own mouthpiece.
- This product is only for oral and dental hygiene. Please do not apply to other areas.
- Do not apply water jet directly under the tongue, inside the ears or insert it into the nostrils. High blood pressure can damage the most sensitive areas.
- In case of failure, anomalies, if the device is dropped, the plug or other parts have impaired or network cable is damaged, do not attempt to repair it yourself. Contact our Customer Service nearest MX ONDA.

CONTROLS AND ELEMENTS

1. Water tank cover
2. Release of the nozzle
3. Water tank or solution
4. Nozzle
5. Switch on / off
6. Pressure selector
7. Charging indicator
8. Micro USB connector for charging the battery

BATTERY (charge and recommendations for use)

Battery dental water jet is shipped with a partial charge and before first use the dental water jet, charge the battery as follows:

Battery charge

1. With the dental water jet in off mode, remove the cover that protects the connector **8** USB Micro for charging the battery.

2. Using the adapter cable (Micro USB - USB) included, connect the dental water jet to a USB port on your computer (which must be on) or use a power adapter with USB output whose voltage is 5 V and is capable of providing 1.5 A (in the case of supplying less current the battery will be slower).
3. During charging, the indicator **7** turn red and when finished charging the battery, the indicator **7** turn blue.
4. Charging is complete, remove the Micro USB connector and replace the cap that protects the connector **8**.

When the battery needs to be recharged?

When you observe pressure decreases drastically or the dental water jet is turned off before two minutes, it will be a sign that the battery is discharged and needs to be recharged as soon as possible.

USE OF DENTAL WATER JET

Filling the water tank

1. With the dental water jet in off mode open the cover **1** of water tank **3**.
2. Fill the tank without exceeding the maximum capacity (which is the transparent tank), you can use weather-water, hot water (without the temperature thereof exceeds 40 °C) or mixed with some kind of mouthwash or mouthrinse liquid.
3. Close the cover **1** of the water tank **3**.

Inserting and removing the nozzle

To insert or remove the nozzle, the dental water jet must be in off mode.

To insert the nozzle, insert it into the hole on the dental water jet by the thickest part of the nozzle and then press it firmly until you hear a "click" indicating that the nozzle is properly retained.

To remove the nozzle press the release button **2** and pull the nozzle.

Using the dental water jet

1. Check that the battery is charged.
2. The dental water jet must be turned off, a full tank of water or solution and the nozzle inserted.
3. Stand over the sink and tilt your head at an angle of about 45 degrees, then enter the nozzle outlet in the mouth.
4. Press the on/off button **5** () and dental water jet will start operating in the last level of the selected flow pressure water.
5. During cleaning, you can increase or decrease the pressure of water flow through the **6** button, each time you press this button will decrease the pressure. If the selected pressure is the minimum, pressing the button **6** again select the maximum pressure.
6. To remove the dental water jet to rinse the mouth, rest or finish cleaning, switch off the irrigator using the on/off button **5** (). Never remove the irrigator of mouth shut if before it.

Very important:

- Do not use this device continuously for more than 5 minutes. Remember that the irrigator has a timer that switches off after 2 minutes.
- To prevent overheating of the pump, never operate the irrigator without water or solution for a long time.

CLEANING

- After each use, empty the water tank.
- Clean the surface irrigator with a damp cloth and dry with another dry. Do not use alcohol or abrasive products.
- Do not insert the device in water or any other liquid.

TECHNICAL SPECIFICATIONS

Water tank capacity	110 ml
Pressure	0.07 ~ 0.55 N
Pulses per minute	1200
Lithium ion battery	3.7 V / 1500 mAh
Charger battery	5 V / 1500 mA max.
Connector for charging	Micro USB
Charging time	4-5 hours approx.
Autonomy	> 2 hours approx.
Dimensions (without nozzle)	64x70x185 mm (L,W,H)
Weight	200 gr

Specifications subject to change without notice.

This product complies with the European Directives 2004/108/EC (EMC – Electromagnetic Compatibility) and the 2006/95/EC (LVD – Low Voltage Directive) relative to Low Tension Security.

“MX, MX ONDA” and its logos are trademarks of MX ONDA, S.A.

Telephone of information and Technical Service: +34 902 551 501

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“MX, MX ONDA” und die entsprechenden logos sind marken der MX ONDA, S.A.